

HTRF® Research Reagents

HTRF AKT1 p-S473 kit - 500 pts

Product Number : 63ADK078PEG Lot Number : 22G

Expiration date : 2026-04-03 Document version : 1

Reagents :	Lot n°
HTRF AKT1 p-S473 kit - 500 pts d2 antibody	22B
HTRF AKT1 p-S473 kit - 500 pts Eu Cryptate antibody	22B
HTRF AKT1 p-S473 kit - Ctrl lysate	25A
HTRF P-T prot. - Lysis Buf.1 (4X) 2 mL	35B
HTRF P-T prot. - Block. reag.(100X) 0.3 mL	28B
HTRF P-T prot. - Block. reag.(100X) 0.3 mL	28B
HTRF P-T prot. - Detect. Buf. 2 mL	32B

Results :

 Incubation : Overnight at RT
 Plate : 384 wells small volume white plate 20 µL
 Reader : Pherastar Plus
 Fitting model : /

Assay characteristics :

	Found value	Specifications Rev 01
Negative	363	/
Control Lysate	3548	/
S/B	9,8	>5
	/	/

HTRF AKT1 p-S473 kit - 500 pts - Lot n° : 22G

 Accepted

 Rejected

Quality control date : 2024-0916	Supervisor's initials : PL
----------------------------------	----------------------------

The information provided in this document is valid for the specified lot number and date of analysis. This information is for reference purposes only and does not constitute a warranty or guarantee of the product's suitability for any specific use. Revvity, Inc., its subsidiaries, and/or affiliates (collectively, "Revvity") do not assume any liability for any errors or damages arising from the use of this document or the product described herein. REVVITY EXPRESSLY DISCLAIMS ALL WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, REGARDLESS OF WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED, ALLEGEDLY ARISING FROM ANY USAGE OF ANY TRADE OR ANY COURSE OF DEALING, IN CONNECTION WITH THE USE OF INFORMATION CONTAINED HEREIN OR THE PRODUCT ITSELF.