

HTRF® Research Reagents

HTRF Human IFN g kit - Std

Product number: 62HIFNGCDA

Lot number: 13A

Product characteristics:

- **Conditioning buffer** Phosphate 50 mM pH7.0 + BSA 1% + stabilizers
- **Appearance** Lyophilized form
- **Storage conditions** 2-8 °C
- **Expiration date** 2028.02.03
- **Reconstitute each vial with** 220 µL distilled H2O

Weight indications:

	Quantity per vial	MW (Da)
Active Moiety (Ab)	4 560.00 pg	N/A

Please note

380µl H2O de volume de reconstitution

The information provided in this document is valid for the specified lot number and date of analysis. This information is for reference purposes only and does not constitute a warranty or guarantee of the product's suitability for any specific use. Revvity, Inc., its subsidiaries, and/or affiliates (collectively, "Revvity") do not assume any liability for any errors or damages arising from the use of this document or the product described herein. REVVITY EXPRESSLY DISCLAIMS ALL WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, REGARDLESS OF WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED, ALLEGEDLY ARISING FROM ANY USAGE OF ANY TRADE OR ANY COURSE OF DEALING, IN CONNECTION WITH THE USE OF INFORMATION CONTAINED HEREIN OR THE PRODUCT ITSELF.