

AlphaLISA®

PPI Buffer (5X)

Product number: AL015C **Lot Number:** 3337899

Product Format: AL015C: 10 mL
AL015F: 100 mL

Manufacturing date: October 7, 2024 **Document version:** 1

Product Information

- Application:** The PPI buffer (5X) has been specially developed for AlphaLISA protein-protein interaction assays.
- Formulation (5X):** PPI Buffer (5X) contains: 250 mM HEPES (pH=7.3), 500 mM NaCl, 2.5% TritonX-100 and 2.5% BSA.
- Storage:** Store PPI Buffer (5X) at 4 °C.
- Stability:** The undiluted stock solutions are stable for at least 6 months from the manufacturing date when stored in original packing at 4 °C.

Please visit our website for additional information on AlphaLISA technology at www.revvity.com

The information provided in this document is valid for the specified lot number and date of analysis. This information is for reference purposes only and does not constitute a warranty or guarantee of the product's suitability for any specific use. Revvity, Inc., its subsidiaries, and/or affiliates (collectively, "Revvity") do not assume any liability for any errors or damages arising from the use of this document or the product described herein. REVVITY EXPRESSLY DISCLAIMS ALL WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, REGARDLESS OF WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED, ALLEGEDLY ARISING FROM ANY USAGE OF ANY TRADE OR ANY COURSE OF DEALING, IN CONNECTION WITH THE USE OF INFORMATION CONTAINED HEREIN OR THE PRODUCT ITSELF.

AL015-R Rev01