

HTRF® Research Reagents

HTRF HP Human IL2 Kit - 500 pts

Product Number : 62HIL02V2PEG Lot Number : 01B


Expiration date : 2026-11-10 Document version : 1

Reagents :	Lot n°
HTRF Diluent 5 (5X) - 12mL	21B
HTRF Detect. Buf. 3 - 1,5 mL	30B
HTRF HP Human IL2 Kit - 500 pts d2 antibody	01A
HTRF HP Human IL2 Kit - 500 pts Eu Cryptate antibody	01A
HTRF HP Human IL2 Kit - Std	01A

Results :

Incubation : 2h at RT
 Plate : 384 wells small volume white plate 20 µL
 Reader : Envision
 Fitting model : 4 Parameter Logistic (4PL 1/y2)* model

h IL2	Conc. (pg/ml)	DF %	DF Accepted Range
Std 1	52	35	24 -- 52
Std 2	126	107	69 -- 121
Std 3	301	251	179 -- 281
Std 4	723	603	428 -- 642
Std 5	1736	1303	984 -- 1475
Std 6	4167	2530	1862 -- 2793
Std 7	10000	2956	2234 -- 3352



Assay characteristics :

	Found value	Specifications Rev 01
Limit Detection	4,2 pg/ml	<52,3pg/ml

HTRF HP Human IL2 Kit - 500 pts - Lot n° : 01B

Accepted
 Rejected

Quality control date : 2024-11-19	Supervisor's initials : PL
-----------------------------------	----------------------------

The information provided in this document is valid for the specified lot number and date of analysis. This information is for reference purposes only and does not constitute a warranty or guarantee of the product's suitability for any specific use. Revvity, Inc., its subsidiaries, and/or affiliates (collectively, "Revvity") do not assume any liability for any errors or damages arising from the use of this document or the product described herein. REVVITY EXPRESSLY DISCLAIMS ALL WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, REGARDLESS OF WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED, ALLEGEDLY ARISING FROM ANY USAGE OF ANY TRADE OR ANY COURSE OF DEALING, IN CONNECTION WITH THE USE OF INFORMATION CONTAINED HEREIN OR THE PRODUCT ITSELF.