

AlphaLISA®

Human Interferon-gamma (IFN- γ) Detection Kit**Product number:** AL217 C **Lot Number:** 3438851**Material provided:****Kit Format:** AL217C: 500 assay points AL217F: 5000 assay points

The number of assay points is based on an assay volume of 50 μ L in 96- or 384-well assay plates using components at the recommended concentrations.

Analyte Format: AL217S: 0.3 μ g (one vial contains an amount of human IFN- γ sufficient for performing 10 AlphaLISA standard curves)**Manufacturing date:** October 10, 2025 **Document version:** 1**Product Information****Application:** This kit is designed for the quantitative determination of human IFN- γ in serum, buffered solution or cell culture medium using a homogenous AlphaLISA assay (no wash steps).**Kit contents:** The kit contains 5 components: AlphaLISA Acceptor beads coated with an Anti-Analyte Antibody, Streptavidin-coated Donor beads, Biotinylated Anti-Analyte Antibody, lyophilized analyte and 10X AlphaLISA HiBlock Buffer.**Sensitivity:** Lower Detection Limit (LDL): 10.6 pg/mL**Dynamic Range:** 10.6 – 100 000 pg/mL**Storage:** Store kit in the dark at +4°C. Store reconstituted analyte at -20°C.**Stability:** This product is stable for at least 12 months from the manufacturing date when stored in its original packaging and the recommended storage conditions. Note: Once reconstituted, the human IFN- γ analyte is stable for at least 60 days at -20°C

Quality Control

Lot to lot consistency is confirmed in an AlphaLISA assay. Maximum and minimum signals, EC₅₀ and LDL were measured on the EnVision Multilabel Plate Reader with Alpha option. We certify that these results meet our quality release criteria. Maximum counts may vary between bead lots and the instrument used, with no impact on LDL measurement.

Maximum signal: 139382 counts

Minimum signal: 765 counts

EC50: 27.17 ng/mL

LDL: 15.049 pg/mL

The information provided in this document is valid for the specified lot number and date of analysis. This information is for reference purposes only and does not constitute a warranty or guarantee of the product's suitability for any specific use. Revvity, Inc., its subsidiaries, and/or affiliates (collectively, "Revvity") do not assume any liability for any errors or damages arising from the use of this document or the product described herein. REVVITY EXPRESSLY DISCLAIMS ALL WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, REGARDLESS OF WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED, ALLEGEDLY ARISING FROM ANY USAGE OF ANY TRADE OR ANY COURSE OF DEALING, IN CONNECTION WITH THE USE OF INFORMATION CONTAINED HEREIN OR THE PRODUCT ITSELF.

AL217-R Rev01

The Revvity logo consists of the word "revvity" in a lowercase, sans-serif font. The letters are black and have a slightly irregular, hand-drawn appearance.

