

Lot specific information

Lot number: 109-25321
 Expiration date: 01-07-2027

Information

The product has been analyzed and is certified to have met the specifications.

Test		Results
Steadylite Plus substrate: 560- (versus Standard Reference): 24505		
Ratio luminescent light-output		0.99
Signal half-life		372 min
Ratio signal half-life		0.96
Background		< 75 cps
Steadylite Plus buffer: 550- (versus Standard Reference): 25141		
Ratio luminescent light-output		0.95
Signal half-life		337 min
Ratio signal half-life		1.04
Background		< 75 cps

This Certificate of Analysis has been produced electronically and is valid without a signature.
 For questions regarding this Certificate of Analysis, please contact your customer service representative.